

# MISCELLANEOUS PRODUCTS



## DOUBLE BUCKET WITH WRINGER



GOOD MANUFACTURING PRACTICES  
CERTIFICATE NUMBER GSCGMP-1-121  
Based on: NOM-251-SSA1-2009

This double bucket allows separating clean water and dirty water\* in different containers making it a very effective tool for highly hygienic cleaning. It has a down press wringer that can be used with wet mops up to 19 oz. that includes an adapter for flat microfiber wet mops (Page 67) and two trays to carry additional cleaning products.

\* Maximum capacity 15 liters each.

PIECE			PACK			PALLET			20 FT CONTAINER		48 FT CONTAINER		53 FT CONTAINER	
CODE	MEASUREMENTS HEIGHT X WIDTH X LENGTH	APPROX. WEIGHT	MEASUREMENTS HEIGHT X WIDTH X LENGTH	CASE QTY.	WEIGHT	CASES	WEIGHT	PIECES	CASES	PIECES	CASES	PIECES	CASES	PIECES
IPC01290	25.60 X 13 X 20 in	17.20 lb	21.26 X 15.75 X 28.15 in	1	20.39 lb	16	360 lb	16	160	160	384	384	416	416

### GENERAL SPECIFICATIONS

SUNLIGHT RESISTANCE		ABRASION RESISTANCE		WATER ABSORPTION		STIFFNESS		ELASTIC RECOVERY	
BUCKET	WRINGER	BUCKET	WRINGER	BUCKET	WRINGER	BUCKET	WRINGER	BUCKET	WRINGER
G	G	A	E	D	D	A	E	E	D

THERMAL RESISTANCE*				CHEMICAL AND MOISTURE RESISTANCE											
DISTORTION		SOFTENING		MOISTURE		DILUTED ACID		CONCENTRATED ACID		DILUTED ALKALI		ALCOHOL AND VEGETABLE OILS		PETROLEUM DISTILLATES	
BUCKET	WRINGER	BUCKET	WRINGER	BUCKET	WRINGER	BUCKET	WRINGER	BUCKET	WRINGER	BUCKET	WRINGER	BUCKET	WRINGER	BUCKET	WRINGER
176 °F	DNA	266 °F	DNA	E	E	E	A	E	D	E	A	E	E	G	G

### SUGGESTED APPLICATION:

For wringing flat or conventional mops, in areas where the risk of contamination is medium-high.

### CONTRAINDICATIONS:

DNA

### SUPPLEMENTARY PRODUCTS:

Olympic One Microfiber frame (IPC00030G), Hygiene Microfiber (IPC01100), Spot Microfiber (IPC01101).



E	EXCELLENT	G	GOOD	A	ADEQUATE (*)	D	DEFICIENT	DNA	DOES NOT APPLY
---	-----------	---	------	---	--------------	---	-----------	-----	----------------

