

SQUEEZES FOR WINDOWS



WINDOW WASHER



GOOD MANUFACTURING PRACTICES
CERTIFICATE NUMBER GSCGMP-1-121
Based on: NOM-251-SSA1-2009

Tool designed for professional window-washing, manufactured with a high-impact polypropylene T handle and a fabric sleeve. Replacement sleeves allow the user to re-use the T handle, resulting in significant long-term savings.

PIECE				PACK			PALLET			20 FT CONTAINER		48 FT CONTAINER		53 FT CONTAINER	
PRODUCT	CODE	MEASUREMENTS HEIGHT X WIDTH X LENGTH	APPROX. WEIGHT	MEASUREMENTS HEIGHT X WIDTH X LENGTH	CASE QTY.	WEIGHT	CASES	WEIGHT	PIECES	CASES	PIECES	CASES	PIECES	CASES	PIECES
WASHER	V04500	18 X 7.09 X 2 in	0.40 lb	6.73 X 10.04 X 18 in	12	5.73 lb	100	487 lb	1200	1000	12000	2400	28800	2600	31200
SLEEVE REPLACEMENT	V85045	18 X 3 X 0.79 in	0.09 lb	5.12 X 9.25 X 10.24 in	10	1.37 lb	285	423 lb	2850	2850	28500	6840	68400	7410	74100

GENERAL SPECIFICATIONS									
SUNLIGHT RESISTANCE		ABRASION RESISTANCE		WATER ABSORPTION		STIFFNESS		ELASTIC RECOVERY	
HANDLE	FABRIC	HANDLE	FABRIC	HANDLE	FABRIC	HANDLE	FABRIC	HANDLE	FABRIC
G	G	G	E	D	E	G	DNA	D	DNA

CHEMICAL AND MOISTURE RESISTANCE											
MOISTURE		DILUTED ACID		CONCENTRATED ACID		DILUTED ALKALI		ALCOHOL AND VEGETABLE OILS		PETROLEUM DISTILLATES	
HANDLE	FABRIC	HANDLE	FABRIC	HANDLE	FABRIC	HANDLE	FABRIC	HANDLE	FABRIC	HANDLE	FABRIC
E	E	E	E	E	G	E	G	E	E	G	D

SUGGESTED APPLICATION:

Ideal for professional window cleaning.

CONTRAINDICATIONS:

DNA

SUPPLEMENTARY PRODUCTS:

Telescopic pole with plastic tip.



E EXCELLENT

G GOOD

A ADEQUATE (*)

D DEFICIENT

DNA DOES NOT APPLY

